

## **v131bd**

→ S2 → DAS profile ba16s\_f.pro → Mars RX (3cm) → Record in parallel on /data/internal → RecordMode "Chans 1,3" → Tape change 317 08:56:00 UT

## **v192a/v192b**

→ Disk (using Huygens cable) → DAS profile av16\_2\_f.pro → Multi RX (21cm) → Record v192a on /data/xraid0 → Record v192b on /data/xraid0 for 3 hours then /data/xraid2 for rest → RecordMode "Chans 1-8" → Start PMDAQ

## **v131be**

→ S2 → DAS profile ba16s\_f.pro → Mars RX (3cm) → Fringe test vc011 → Record in parallel on /data/internal → RecordMode "Chans 1,3" → Tape change 318 21:56:00 UT → Tape change 319 03:56:00 UT

## **Disk change to Xraid W1-14**

## **vt11b1**

→ Who knows TBA by Chris

## **Disk change to Xraid X1-7**

## **v190e**

→ Disk (using Huygens cable) → DAS profile bv16\_2\_f.pro → Mars RX (3cm) → Fringe test vc013 (Record on /data/removable) → RecordMode "Chans 1-8" → Record on /data/xraid0 for 6 hours → Record on /data/xraid1 for 6 hours → Disk change ~15 UT

## **v215a**

→ Disk (using normal cable) → DAS profile aa4s\_n.pro → Methanol RX (5cm) → Fringe test vc014 (Record on /data/internal) → RecordMode "Chans 1,3" → Record on /data/removable

## v131bf

→ S2 → DAS profile ba16s\_f.pro → Mars RX (3cm) → Fringe test vc015 → Record in parallel on /data/internal → RecordMode “Chans 1,3” → Tape change 321 21:56:00

## v190d

→ Disk (using Huygens cable) → DAS profile av16\_2\_f.pro → H-OH RX (18cm) → Fringe test vc016 (Record on /data/removable) → RecordMode “Chans 1-8” → Record on /data/xraid0 for 6 hours → 15:30 UT disk change → Record on /data/xraid1 for 6 hours → 21:30 UT disk change → Record on /data/xraid2 for 6 hours → 03:30 UT disk change → Record on /data/xraid3 for 6 hours

From:

<https://www.atnf.csiro.au/people/Tasso.Tzioumis/vlbi/dokuwiki/> - **ATNF VLBI Wiki**

Permanent link:

<https://www.atnf.csiro.au/people/Tasso.Tzioumis/vlbi/dokuwiki/doku.php/lbaops/parkes/nov2006>

Last update: **2015/12/18 16:39**

